

ZOLL AED 3® PADS

CPR Uni-padz® III Universal (Adult & Pediatric)

REV071624

Truly universal, these pads can be used in an adult or child rescue. For child rescue, reposition the CPR Uni-padz as shown and activate the child rescue mode. CPR Uni-padz electrodes provide full CPR feedback for adult patients and have a five-year installed shelf life.


Every set of CPR Uni-padz electrodes comes with a Rescue Accessory Kit that includes scissors for cutting clothing to expose the victim's chest, non-latex gloves and a breathing mask for rescuer protection, a paper towel to ensure dry contact with electrodes, a razor for removing hair, and a wipe for cleanup.

Features:

- **INCLUDES MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE**
- Compatible with the ZOLL® AED 3 Automatic Defibrillator
- For use with adult patients and (in child rescue mode) for children under the age of 8 years or less than 55 pounds (25 kilograms)
- Provides full CPR feedback for adult patients
- Has a five-year shelf life

Dimensions:

- Package: L 8.5 in. x W 8.24 in. x D 1.25 in.
- Packaged Weight: 12.4 oz.



ZOLL CPR Uni-padz® III

ITEM#	NSN#
75-0027	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

PRODUCTS
WITH A **MISSION®**