

ZOLL AED PLUS® PADS

Pedi-padz® II Electrodes (One Pair)

REV071624

ZOLL's intelligent pediatric capability is built into the AED Plus. The AED Plus Defibrillator (sold separately) can sense when the Pedi-padz II Electrodes are being used, and will switch the AED Plus from adult to pediatric mode, modifying the device so it will provide an appropriate pediatric ECG analysis algorithm and therapy energy level. This product is for use with patients under the age of 8 years or less than 55 pounds (25 kilograms)* and has a two-year shelf-life.


**Caution: Pediatric electrodes should never be used with an adult patient*

Features:

- **INCLUDES MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE**
- Compatible with the ZOLL® AED Plus Automatic Defibrillator
- For use with patients under the age of 8 years or less than 55 pounds (25 kilograms)
- Has a two-year shelf life

Dimensions:

- Package: L 8.5 in. x W 8.25 in. x D 1 in.
- Packaged Weight: 7 oz.



ZOLL PEDI-PADZ® II

ITEM#	NSN#
75-0022	----



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

PRODUCTS
WITH A **MISSION®**