

ENT POCKET LIGHT SET

REV060524



Lightweight and Durable ENT Diagnostic Sets for Field Care...

The ENT Pocket Light Set is designed to provide bright, white illumination for effective cerumen management. Its versatility, durability, and ease of use make it ideal for use in both clinical and field settings.

This set includes a curette for removal of debris from the ears, eyes, and nose. The light bender is used to provide easy access to the nasal cavity.

Features

- Compactable for field use
- Made for easy cerumen removals with our flexible nylon for ears, nose and eyes
- Bright white illumination from our pocket light for excellent field use scenarios
- Added capabilities for exploring the nasal cavity and ears with the light bender and speculum

Specifications

- Packaged Dimensions: L 8.5 in. x W 5 in. x D 1 in.
- Unpacked Dimensions: L 6.75 in. x W 3.25 in. x D 1 in.
- Weight: 0.6 lbs

ENT Pocket Light Set

ITEM#	NSN#
20-0118	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**