

# NAR-5 SAR AID KIT™

**The standard casualty response bag which meets the requirements for the SAR Level A Medical Bag**

The NAR 5 Aid Kit™ from North American Rescue® is designed to meet the rigorous requirements of the U. S. Navy Search and Rescue (SAR) Swimmers. This ergo dynamic, low profile & durable casualty response bag is constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self repairing zippers with silent 550 cord pulls and is **designed to be hung in aircraft or other transport vehicle by top bag handle**. The NAR 5 also features 6 hook-and-loop internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products.

Pouches are designed for quick access to medical equipment necessary to respond to specific life threatening medical emergencies. Pouches attach via hook-and-loop to allow the first responder to remove the specific pouch needed to provide on-scene medical care to the patient.

Slim-line, low-profile bag is securable using anti-pilferage seals to protect all medical components inside the bag. Kit can be carried as a backpack using shoulder straps that stow away until needed.

## FEATURES/SPECIFICATIONS

- Meets the requirements for the SAR Level A Medical Bag
- Ergo dynamic, low profile & durable casualty response bag constructed of Berry Amendment 500D CORDURA® solution dyed nylon
- GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, selfrepairing zippers with silent 550 cord pulls
- Featuring 6 hook-and-loop internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products

## DIMENSIONS

- Packaged: H 24 in. x W 15 in. x D 7 in.
- Weight: 14 lbs. 10 oz.




Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)

## NAR-5 AID KIT <sup>Rx</sup> ONLY

ITEM #	OPTION	NSN#
80-0416	Kit	---
80-0399	Bag Only	8465-01-669-5796*

\*NATO NSN

# NAR-5 SAR AID KIT™



## KIT CONTENTS

- 1 x Bag (CCRK® - NAR-5 SAR Aid)
- 4 x NAR® PPE Kit™ (Personal Protection Equipment Kit)
- 4 x C-A-T® (Combat Application Tourniquet®)
- 4 x Combat Gauze® Z-fold Hemostatic
- 4 x NAR® Compressed Gauze
- 4 x ETD® 4 in. Emergency Trauma Dressing
- 1 x ETD® Abdominal Emergency Trauma Dressing
- 4 x NAR® Nasopharyngeal Airway - Pre-lubricated, 28F
- 4 x ARS® for Needle Decompression (14G x 3.25 in.)
- 4 x HyFin® Vent Twin Pack Chest Seal
- 1 x Cyclone® New & Improved BVM
- 1 x CPR Mask
- 4 x Triangular Bandage
- 4 x 3 in. Elastic Wrap Bandage
- 2 x Cervical Collar
- 2 x SAM® Splint II
- 1 x SAM® Pelvic Sling
- 1 x Roll Surgical Tape (1 in.)
- 1 x NAR® Trauma Shears (7.25 in.)
- 1 x Sterile Burn Sheet (60 x 96 in.)
- 12 x Flexible Fabric Bandages
- 4 x Polycarbonate Eye Shield (PES)
- 1 x NAR® Survival Blanket
- 2 x 1000cc IV Solution
- 1 x Mini Maglight® with batteries

*Designed to meet the rigorous requirements of the U. S. Navy Search and Rescue (SAR) Swimmers*



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277