

FIELD CORPSMAN/ENT ACCESSORIES

REV062524

Lightweight and Durable ENT Diagnostic Accessories for Field Care...

Replacement accessories/components compatible with the Basic and Deluxe Field Corpsman kits. These components fit over the head of the ENT Pocket Exam Light.

Dimensions

Packaged: 4 in. x 6 in.



ENT Otoscopic Head

Wt:
0.25 lb.

ITEM#	NSN#
ZZ-1674	----

- Used to gain better visibility into the ear canal



Speculum Adaptor (For otoscopic head)

Wt:
0.2 lb.

ITEM#	NSN#
ZZ-1675	----

- Fits onto ENT Otoscopic Head to attach disposable specula



Nasal/Head Light Bender

Wt:
0.3 lb.

ITEM#	NSN#
ZZ-0042	6515-01-491-2101

- Used to look in nose, test pupillary reflex, or focus light into the ear or mouth



Nylon Eye Loop (10pk)

Wt:
0.75 lb.

ITEM#	NSN#
20-0012	6515-01-491-2147

- Designed to gently sweep foreign bodies from the surface of the eye.
- Lighted tip illuminates work area
- Comes in a pack of 10



Lighted Stainless Steel Curette (10pk)

Wt:
0.75 lb.

ITEM#	NSN#
20-0020	6650-01-491-2008

- Used to help remove wax from ear, or assist in debris removal from small wounds
- Lighted tip illuminates work area
- Comes in a pack of 10



Eye Magnet (10pk)

Wt:
1.14 lb.

ITEM#	NSN#
20-0022	6545-01-491-0855

- Helps to remove debris from inside the eyes
- Lighted tip illuminates work area
- Comes in a pack of 10



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**