

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT *Nylon*

REV062824

Control the bleed with these easy-to-use tools!

Life threatening bleeding can occur from everyday emergencies such as vehicle crashes, incidents involving farming equipment or public transportation, accidental injuries from natural disasters like tornadoes, and inflicted injuries from man-made disasters such as shootings and bombings. The Individual Bleeding Control Kit by North American Rescue® is designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding.

Contained in a nylon carry bag, this life saving kit provides intuitive and easy-to-use tools that are proven to help to save lives such as the fastest, safest and most effective prehospital tourniquet in the world, the C-A-T®, quality gauze, shears, and pressure dressings/bandages. Included are step-by-step illustrated “just in time” pictorial instructions.

NAR® Individual Bleeding Control Kits are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced CELOX™ configurations.

BASIC:

The NAR Individual Bleeding Control Kit provides life saving bleeding control equipment such as tourniquets, pressure dressings and gauze bandages and puts them into the hands of the public to help control bleeding and save lives.

INTERMEDIATE:

Includes all equipment found in the Basic Individual Bleeding Control Kit with the added capability to treat penetrating chest wounds with the HyFin® Vent Compact Chest Seal Twin Pack (for entrance and exit wounds).

ADVANCED / ADVANCED BCD / ADVANCED CELOX™:

Includes all equipment found in the Intermediate Individual Bleeding Control Kit with the added capability to treat major bleeding with a hemostatic dressing that enhances blood clotting abilities—for faster bleeding control.

Features:

- Designed to provide essential equipment that empowers the general public to take action as immediate responders in stopping life threatening bleeding
- Provides intuitive and easy-to-use tools that are proven to help to save lives
- Kits are available in Basic, Intermediate, Advanced, Advanced BCD, and Advanced Celox™ configurations
- Packaged in a compact, red nylon carry bag

Specifications:

- H 7 in. x W 5 in. x D 2.5 in.
- Weights:
 - Basic: 15 oz
 - Intermediate: 16.5 oz
 - Advanced: 17.4 oz
 - Bag Only: 3.41 oz



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT *Nylon*

REV062824

Every PUBLIC ACCESS INDIVIDUAL BLEEDING CONTROL KIT

includes these core items*:



QTY 1 C-A-T® Tourniquet

QTY 1 Responder 6 in. ETD

QTY 2 NAR® Gauze, Wound Packing

PAIR 2 Responder Nitrile Gloves, Lg

QTY 1 Permanent Marker, Small

QTY 1 Responder Trauma Shears, Lg

QTY 1 Instructions Bleeding Control

QTY 1 NAR® Survival Blanket

Individual Bleeding Control Kit

Level	Item#
Basic	80-0452
Intermediate	80-0453
Advanced	80-0454
Advanced with BCD	80-0523
Advanced with CELOX™	80-1570
Bag Only	80-0459W

PLUS these Individual Components based upon your configuration of choice:



QTY 1 HyFin® Vent Compact, Twin



QTY 1 HyFin® Vent Compact, Twin



QTY 1 Combat Gauze® LE



QTY 1 HyFin® Vent Compact, Twin



QTY 1 Bleeding Control Dressing



QTY 1 HyFin® Vent Compact, Twin



QTY 1 5' CELOX™ Hemostatic

All these kits include:

- 1 x C-A-T® Tourniquet
- 1 x 6 in. NAR® Responder ETD® (Emergency Trauma Dressing)
- 2 x NAR® Wound Packing Gauze
**Advanced Kits replace (1) Wound Packing Gauze with (1) Hemostatic Gauze.*
- 2 x Pair, Responder Nitrile Gloves, Lg
- 1 x Responder Trauma Shears, 7.25 in.
- 1 x NAR® Survival Blanket
- 1 x Permanent Marker, Small
- 1 x Just in Time Instruction Card
- 1 x PVC Bleeding Control Patch
- 1 x Red Nylon Bag

Intermediate Kit also includes:

- 1 x HyFin® Vent Compact, Chest Seal - Twin Pack

Advanced Kit also includes:

- 1 x HyFin® Vent Compact, Chest Seal - Twin Pack
- 1 x Combat Gauze® LE, Z-Fold Hemostatic Bandage

Advanced with BCD Kit also includes:

- 1 x HyFin® Vent Compact, Chest Seal - Twin Pack
- 1 x QuikClot® Bleeding Control Dressing

Advanced with CELOX Kit also includes:

- 1 x HyFin® Vent Compact, Chest Seal - Twin Pack
- 1 x CELOX™ Z-Fold Hemostatic Gauze



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**