




# NAR® NASOPHARYNGEAL AIRWAY

28F (117 MM) - PRE-LUBRICATED

REV071124



**NAR® Nasopharyngeal Airway, 28F - Pre-lubricated**

QTY	ITEM#	NSN#
Single	10-0060	6515-01-710-9764
Box of 10	10-0067	6515-01-710-9783

## Introducing the North American Rescue® Pre-Lubricated Nasopharyngeal Airway

Introducing the North American Rescue Pre-Lubricated Nasopharyngeal Airway

Our **Made in USA** pre-lubricated nasopharyngeal airway is an airway adjunct that can be used to secure an airway on patients with or without an intact gag reflex. The beveled tip and pre-lubricated, non-latex, soft flexible material ensures safe and effective insertion. Easy to use and ready to insert after opening, the pre-lubricated 28F (117mm) nasal airway provides maximum comfort and ease of insertion. Packaged in a durable pouch, the pre-lubrication provides smooth insertion, minimizing discomfort and potential trauma while eliminating procedure prep time. Additionally, the anatomical shape and soft material offer optimal conformity, reducing the risk of tissue damage and enhancing patient safety. Made in USA with high grade materials and designed for optimal efficiency and ease of use.

### Specifications:

- Packaged pre-lubricated to decrease procedure prep time
- Water based lubricant (-50°F freezing point)
- Low cube and weight
- Rounded bevel tip allows for gentle, optimal insertion
- Soft material avoids damage to nasal passage
- Unique proprietary material and design maintains patent shape
- Durable packaging
- Sterile
- Latex free
- Made in USA

### Dimensions:

- Ind. Package: L 7.5 in. x H 3.5 in.
- Box: L 9 in. x W 4 in. x D 4 in.
- Device Size: 28 FR (117 mm)
- Weight: 0.02 lbs (Ind.)  
0.50 lbs (Box)



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**