

# NAR® TRAUMA SHEARS

REV071224



6.25 in. Small

7.25 in. Large  
featuring O<sub>2</sub> wrench

*Unique, all-over black  
Trauma Shears ideal for the  
tactical environment*

## Premium-grade stainless steel with precision sharpness...

Built for safety and efficiency, the NAR Trauma Shears from North American Rescue let tactical healthcare professionals get the job done, and get it done right.

These paramedic shears are made from premium-grade stainless steel. Their serrated lower edge cuts through the toughest materials. Their soft-grip, non-slip contoured handles provide maximum control, even when wet. Their clothing lift lip and rounded tip provide the ultimate in patient safety.

These are black shears made to meet the challenges of casualty care with precision and dependability. Our medical scissors come in two convenient lengths, 6.25 or 7.25 inches.

### Features:

- Premium-grade stainless steel with precision sharpness
- Enhanced serrated lower edge cuts through the toughest materials
- Clothing-lift lip and rounded tip for ultimate patient safety
- Soft-grip, non-slip contoured handles for maximum control even when wet
- Unique all-over black color
- Large shears include an O<sub>2</sub> wrench

### Dimensions:

- Small:  
L 6.25 in. x W 2.75 in. x D 0.3 in.  
Weight: 1.5 oz
- Large:  
L 7.25 in. x W 3.7 in. x D 0.4 in.  
Weight: 2 oz

### NAR® Trauma Shears

| SIZE  | ITEM#   | NSN#              |
|-------|---------|-------------------|
| Small | ZZ-0064 | 6510-01-622-1624* |
| Large | ZZ-0063 | 6515-01-538-9276  |

\* NATO NSN#



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**