

CRÉDO PROMED™ SERIES FOUR

472 THERMAL INSULATED CHAMBER

REV071224



Proven Payload Protection...

The Crêdo ProMed™ Series 4 - 472 Thermal Insulated Chambers are inner cooling panels that are designed to fit inside of the Crêdo ProMed™ Series Four Containers, which are sold separately. These products work together to thermally protect the integrity of valuable medical supply payloads for 72 - 96 hours.

The patented Thermal Insulation Chamber (TIC™) panels are constructed from phase change material and VIP components contained within, and are qualified to consistently protect medical materials such as blood and platelets within 2-8°C.

Easy conditioning of the TIC™ system panels makes the ProMed a simple, nimble, and consistently reliable option. Once the TIC™ system panels are staged and properly placed inside the Vacuum Insulated Panel (VIP) assembly of the Crêdo ProMed™ Series Four Container (sold separately), simply insert payload content and be on your way.

Please note that the 2 Liter Thermal Insulated Chamber is a container with a separate lid, while the 4 Liter and 8 Liter version consists of separate panels that can be configured into a cube shape.

Features:

- Duration: 48 hours
- Temperature: 2° to 8° C
- Volume: 2L, 4L, 8L

Dimensions:

2 Liter

- 7.5 in. L x 6.5 in. W x 5.75 in. D
- Weight: 3 lbs 12.2 oz

4 Liter

- Per panel: 6 in. x 6 in. (6 pcs)
- Weight: 5 lbs 3.2 oz

8 Liter

- Shorter Panels: 6.5 in. x 6.5 in. (2 pcs)
- Longer Panels: 6.5 in. x 11 in. (4 pcs)
- Weight: 7 lbs 1.6 oz

Crêdo ProMed™ Series Four 472 Thermal Insulated Chamber

CAPACITY	ITEM #
2L	35-0021
4L	35-0022
8L	35-0023



Container
Sold
Separately



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**