

M-FAK LCL Bag Only

REV072324

NAR's most compact IFAK bag, now available in laser cut nylon!

The Mini First Aid Kit (M-FAK) bag - now featuring Laser Cut Laminate (LCL) MOLLE - is a super compact, rugged nylon bag that opens via 550 cord zipper pulls from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 inch internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet (sold separately). Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips on the front, and (2) two MOLLE strips on the spine's edge to attach additional gear such as TQ/holders, tactical lights, or knives (sold separately) if needed.

Features:

- NAR's most compact, versatile Mini First Aid Kit bag - now with laser cut laminated MOLLE
- The LCL IFAKS are double spaced for the MOLLE attachment points on the rear of each IFAK for maximum adjustability and super secure mounting to any platform
- Half the size of our popular T.O.R.K. platform

Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 8 oz

M-FAK LCL BAG

COLOR	ITEM #
BLK	80-1601
COY	80-1603
OCP	80-1604
ODG	80-1670
RGR	80-1602
RED	80-1610



Nylon Colors:

BLK

COY

OCP

ODG

RGR

RED



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**