

# APRU ACCESSORIES

AUTONOMOUS PORTABLE REFRIGERATION UNIT

REV082224

## APRU Custom Foam Transport Case

The APRU Custom Foam Transport Case is a heavy duty plastic transport box with laser-cut, custom foam for optimal protection of the Autonomous Portable Refrigeration Unit (APRU) and its accessories.

### Features:

- Heavy duty, with laser-cut custom foam inserts
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) and its accessories
- APRU and its accessories are sold separately

### Dimensions:


- Dimensions  
- H 30 in. x W 18 in. x D 23 in.
- Weight: 41 lbs.



*Contents  
Sold  
Separately*

### Transport Case

ITEM #	NSN #
35-0035	----



### Vehicle Mount & Cable

ITEM #	NSN #
35-0036	----

## APRU Vehicle Mount & Hardware Cable

The APRU Vehicle Mount and Hardwire Cable stabilizes the Autonomous Portable Refrigeration Unit (APRU) when transported in a vehicle. Includes hardwire cable and installation support. Does not include strap.

### Features:

- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately
- Does not include strap

### Dimensions:

- Vehicle Mount Dimensions:  
H 1.25 in. x W 10 in. x L 19.25 in.
- Hardwire Cable Length: 15 ft.
- Total Weight: 5.2 lbs.



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# APRU ACCESSORIES

AUTONOMOUS PORTABLE REFRIGERATION UNIT

## APRU Backpack Carrier with Pouches

The APRU Backpack Carrier with Pouches is the CRO Medical (CRO) carrier for the Autonomous Portable Refrigeration Unit (APRU), adding carrying capability by hand or as a backpack. Also provides additional impact protection. Includes 2 extra battery pouches and a power cable storage pouch.

### Features:

- Adds carrying capability by hand or as a backpack
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately
- Includes two extra battery pouches and a power cable storage pouch

### Dimensions:

- Dimensions: H 19 in. x W 9.5 in. x L 10.5 in.
- Weight: 3.2 lbs.



### Backpack Carrier

ITEM #	NSN #
35-0037	4110-01-706-0999



REV082224

### Soft Shoulder Carrier

ITEM #	NSN #
35-0038	----

## APRU Soft-Shoulder Carrier (FTECH)

The APRU Soft-Shoulder Carrier is compatible with the Autonomous Portable Refrigeration Unit (APRU), adding carrying capability by hand. Also provides additional impact protection.

### Features:

- Adds carrying capability by hand
- Provides additional impact protection
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately

### Dimensions:

- Dimensions: H 19.5 in. x W 10 in. x L 11 in.
- Weight: 3 lbs.



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**