

DELTA ICE TEMPERATURE CONTROL INSERT 2 LITER

REV080224

Optimal temperature storage in virtually any environment.

The Delta ICE's Temperature Controlled Insert is an aluminum insert featuring integrated, built-in NIST traceable sensors that promote temperature uniformity and an even thaw, eliminating hot spots. This SMART technology maintains two (2) 500ml of whole blood for 48 hours at 2-6°C and eliminates the need for glycol.



Features:

- For use with the Delta ICE Smart Blood Cooler (sold separately)
- Aluminum insert featuring integrated, built-in NIST traceable sensors
- Maintains two 500 ml of whole blood for 48 hours at 35.6 - 42.8°F, and eliminates the need for glycol

Dimensions:

- Dimensions: H 5.75 in. x W 6.25 in. x D 8.24 in.
- Weight: 5 lbs.

DELTA ICE TEMPERATURE CONTROL INSERT - 2L

ITEM #	NSN #
35-0043	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**