

SMART BLOOD COOLER

DELTAICE_{2L}



« Explore this product
at NARescue.com



REMOTE ALERTS

Bluetooth and Cellular alerts warn you and your command of temperature excursions

FDA LISTED

Demonstrate safety simply by showing our FDA listing to your blood supplier

RUGGED DESIGN

Insulation protected from drops and punctures

IN REAL-TIME

Real-time temperature reading and recording

SMART INTEGRATION

Temperature-controlled, aluminum insert with intergrated traceable sensors promote temperature uniformity and ease of use

KNOW YOUR STORAGE TIME

Screen and app display remaining time before thermal storage melts

RECHARGEABLE BATTERY

Recharge with USB type C

LONG STORAGE TIME

Stays cold up to 48 hours, while internal memory retains data for 3 days

EXTREME HEAT ENVIRONMENTS

Stays cold 12 hours even in uncooled vehicles



REV071824



DISTRIBUTED BY

NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

SMART BLOOD COOLER DELTAICE 2L

SPECIFICATIONS



GENERAL

Size

11.5"L x 9.7"W x 9.5"H

Weight

10.5 lbs unloaded

FDA Classification

Class II. Code KSE

Standard

Meets AABB Standards for Blood Banks and Transfusion Services

STORAGE

Duration

48 hours in 25°C (77°F)

24 hours in 44°C (111.2°F)

12 hours in 60°C (140°F)

Internal Volume

2 L

Unit capacity

Two 500 ml bags



ENVIRONMENTAL

Temperature Ranges:

Operating

2° to 60°C (36°F to 140°F)

Storage

-20° to 70°C (-4°F to 158°F)

Ingress Protection

IP55

TEMPERATURE RECORDING

Frequency

15 minutes

Duration

3 days

Accuracy

± 1°C

Reporting

Cellular and Bluetooth to Application

PREPARATION

Freezer Space

4.7 L

Freeze Time

8 hrs at -9.5°C (15°F)

Prep Time

15 minutes in 25°C (77°F) ambient

BLACK: REF 35-0040

BLUE: REF 35-0041

RED: REF 30-0042



« Explore this product at NARescue.com



DISTRIBUTED BY

NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

35 Tedwall Ct. • Green, SC 29651 • USA
Info@NARescue.com

DELTA
DEVELOPMENT TEAM®