

# WIZARD WAND EYE CARE KIT

REV082224

## Lightweight and Durable ENT Diagnostic Sets for Field Care...

The Wizard Wand Eye Kit is an innovative eye care product designed for the examination and immediate care of corneal abrasions and foreign bodies within the eye. The Wizard Wand's cobalt blue light source is a collimated, diffused, narrow spot of light which prevents shadows and irregularities in the illuminated area. The cobalt blue light effectively excites standard fluorescein dyes (Fluorescein strips not included) in order to locate corneal abrasions and translucent particles on the surface of the eye.

Two attachments are included with the Wizard Wand. The lighted nylon eye loop facilitates the removal of non-metallic particles from the surface of the eye and the conjunctiva. The lighted eye magnet effectively removes non-embedded metallic particles without contacting the surface of the eye. The Wizard Wand is powered by two AAA batteries.

### Features

- Eye care management set
- Compact eye care kit for field use
- Built to remove foreign bodies in the eye
- Fluorescein strips not included

### Specifications

- Packaged Dimensions: L 8.5 in. x W 5 in. x D 1 in.
- Unpacked Dimensions: L 6.75 in. x W 3.25 in. x D 1 in.
- Weight: 0.6 lbs.

### Wizard Wand Eye Care Kit

ITEM#	NSN#
20-0100	----



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A **MISSION®**