

BREAKAWAY MEDICAL SIMULATION UNIFORM SET

REV082324

The Breakaway Medical Simulation Uniform Set comes with one shirt and one pair of pants, and comes in Scorpion OCP (Operational Camouflage Pattern). The set is one size - XL.

The Breakaway Medical Simulation Pants features side cargo pockets, and multiple hook and loop fasteners, from waist to bottom of crotch and down the crotch inseam of both legs.

The Breakaway Medical Simulation Shirt features two breast pockets, with loop piling name tags over the pockets. There is a hook and loop closure down center of chest, closures down the length of the sleeves, and down the sides from armpit to waist. All of the hook and loop closures allow the garment to breakaway without cutting for repeated use.

Sized to fit both TOManikin and Tamikin simulators (not included), the hook and loop closures simplify dressing the manikin and also allow trainees the ability to quickly remove garments during patient assessment scenarios.



Specifications:

- Breakaway Medical Simulation Uniform Set includes one shirt and one pair of pants
- Features multiple hook and loop closures, allowing the garment to breakaway without cutting for repeated use
- One size - XL
- Fits TOMManikin and Tamikin Simulators (not included)

Breakaway Medical Simulation Uniform Set

ITEM#	OPTION	NSN#
93-1018	Set	----
93-1019	Shirt	----
93-1020	Pants	----

Dimensions:

- *Packaged Measurements:*
L 16.5 in. x W 12 in. x D 3 in.
- *Weight:* 3 lbs.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

BREAKAWAY MEDICAL SIMULATION UNIFORM SET

REV082324



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**