

# ROO® M-FAK LCL *Bag Only*

REV092424

**One of our most compact bags, now available in laser cut nylon!**

The ROO M-FAK's super compact nylon bag features a unique front tourniquet pouch that easily opens with the pull of a tab for rapid one handed access to your tourniquet (sold separately) without the need to open the full kit. Multiple elastic loops on the interior allow for secure gear organization, while it's rugged exterior includes a Laser Cut Laminated (LCL) MOLLE style attachment for mounting in the standard vertical position as well as a 3 in. internal sleeve that accommodates horizontal attachment to your duty belt.

## Features:

- Bag ONLY, contents not included
- One of NAR's most compact, versatile Mini First Aid bags featuring rapid access front pouch
- 500 Denier nylon bag with interior elastic loops for secure gear organization
- Super compact, rugged nylon platform that allows attachment both vertically and horizontally
- The LCL ROO IFAKS bags are double spaced for the MOLLE attachment points on the rear of each IFAK bag for maximum adjustability and super secure mounting to any platform
- Compatible with NAR® Quick Detach Panels (sold separately)

## Specifications:

- L 7 in. x W 3.75 in. x D 4 in.
- Weight: 8 oz

## ROO® M-FAK LCL Bag

COLOR	ITEM #
BLK	80-1605
COY	80-1607
OCP	80-1608
ODG	80-1646
RED	80-1611
RGR	80-1606



TOURNIQUET  
NOT  
INCLUDED

Nylon  
Colors:

BLK

COY

OCP

ODG

RED

RGR



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**