

TRAUMA & FIRST AID KITS *RESUPPLY - Class A*

REV082924

The workplace First Aid Kit that exceeds ANSI/ISEA Z308.1-2021 for Class A

The Trauma and First Aid Resupply Kit - Class A from North American Rescue® is designed as a portable First Aid Kit for the workplace that exceeds ANSI/ISEA Z308.1-2021 Standard for Class A kits. Packaged in a convenient zippered nylon bag that is easily transported to the patient's side, this kit offers the ability to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T® Tourniquet, the safest and most effective tourniquet on the market today. Additionally, kit options are available with one of two specialized hemostatic gauzes: QuickClot® Bleeding Control Dressing, or Celox™ First Aid Hemostatic Dressing, that helps rapidly control bleeding.

Features:

- Workplace first aid kits that exceed ANSI/ISEA Z308.1-2021 Standard for Class A kits
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the safest and most effective Tourniquet on the market
- Designed so the kit can be easily carried to the patient's side

Dimensions:

- L 15.0 in. x W 12.0 in. x D 2.75 in.
- Weight: 2.35 lb



Trauma & First Aid Resupply Kits - Class A

ITEM #	OPTION
80-1153	w/ Compressed Gauze
80-1155	w/ Bleeding Control Dressing
80-1479	w/ CELOX™ Gauze




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

TRAUMA & FIRST AID KITS *RESUPPLY - Class A*

REV082924




Kit Contents:

- 1 x C-A-T® Tourniquet, Orange
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x NAR® Responder Compressed Gauze
- 1 x 4 in. NAR® Responder ETD®
- 1 x 4 in. Elastic Wrap Bandage
- 1 x 2 in. Elastic Wrap Bandage
- 1 x Triangular Bandage
- 1 x SAM® Splint
- 2 x Gauze, Pad - 4 in. x 4 in. (Pack of 2)
- 2 x Gauze Pad - 5 in. x 9 in.
- 16 x Flexible Fabric Bandages, 1 in. x 3 in.
- 1 x BurnTec® Dressing - 4 in. x 4 in.
- 1 x CPR Microshield Mask
- 10 x Antiseptic Towelette
- 10 x Antiseptic Ointment
- 1 x Roll 1 in. Surgical Tape
- 1 x Instant Cold Pack - 6 in. x 9 in.
- 1 x Eye Wash Solution, 1 oz
- 2 x Sterile Oval Eye Pad
- 2 pr Responder Nitrile Gloves - Large
- 1 pr Responder Trauma Shears - Large
- 10 x Hand Sanitizer
- 10 x Burn Gel
- 1 x NAR® Survival Blanket
- 1 x Black Permanent Marker, Large
- 1 x First Aid Pocket Guide
- 1 x Just In Time Bleeding Control Instructions
- 1 x Trauma & First Aid Kit Contents Inventory Card
- 1 x Inspection Card

Plus one of these:

- 1 x NAR® Responder Compressed Gauze
- 1 x QuikClot® Bleeding Control Dressing
- 1 x Celox™ First Aid Hemostatic Dressing

PLUS+
one of these:



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**