

MW EVA LITTER 23S *With Case*

REV090624



Specifications:

- Folded: H 17 in. x W 7.5 in. x D 5 in.
- Open: H 72.8 in. x W 20 in. x D 2 in.
- Weight of litter (poles and bed): 6.37 lbs.
- Weight (with bag): 7.4 lbs.
- Max carrying load: 330.7 lbs. (165 lbs. per pole)

For Particularly Demanding Environments

The MediWear EVA Litter 23S is a quick, easy to assemble lightweight, compact casualty litter that allows rapid transport of casualties. The telescopic poles are constructed out of extremely durable carbon fiber, with new aluminum pole ends that quickly and easily deploy. The litter includes a carrying case that can be secured to load bearing equipment using the MOLLE attachment system.

The light weight and small size of the product make it suitable for carrying into particularly demanding environments where carrying capacity, muscle strength, and reaction time can be the difference between a life lost or a life saved.

Features:

- Weighs 6.37 lbs. and can carry up to 330.7 lbs. (165 lbs. per pole)
- Telescopic poles are constructed of carbon fiber with new aluminum carry handles at each pole end, making it significantly easier to deploy
- Has integrated patient belt
- Litter Bed consists of 75% PA and 25% PU
- Fabric finish is water and dirt repellent

MW EVA LITTER 23S

ITEM#	NSN#
60-0093	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

MW EVA LITTER 23S *With Case*

REV090624



**NAR-4 Aid Bag
sold separately**



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**