

TORK® LCL (Tactical Operator Response Kit)

REV112624



Point-of-wounding IFAK for the Assault Vest or Body Armor, now available in laser cut vinyl!

The Tactical Operator Response Kit (TORK) by North American Rescue® - *now featuring Laser Cut Laminate (LCL) MOLLE* - is designed to meet the need for an individual first aid kit carrying basic life support equipment, in order to address leading causes of potentially preventable death in a tactical environment.

This TORK® contains essential life-saving equipment for point-of-wounding, self-aid/buddy-aid to address injuries caused by penetrating or blast trauma. Its sleek, streamlined size allows for easy attachment using MOLLE/PALS style connectors to the assault vest/body armor or NAR's Quick Detach Panel (available for separate purchase). The TORK® also features a signature Red-Tip Technology® rapid opener for fast and easy, one-pull tab opening of the kit.

The TORK® LCL is available in different levels of care: Basic, Basic with Combat Gauze, and Advanced with Combat Gauze.

Features:

- Essential life-saving tools to address injuries caused by penetrating or blast trauma
- Streamlined Individual First Aid Kit (IFAK) providing Basic Life Support for immediate Self-Aid/Buddy-Aid at the point of wounding
- Includes the added advantage of a hemostatic dressing for extra capability in the control of life threatening bleeding
- The LCL IFAKS are double spaced for the MOLLE attachment points on the rear of each IFAK for maximum adjustability and super secure mounting to any platform
- Also attaches to NAR's Quick Detach Panel (sold separately)

Specifications:

- Dimensions: H 7 in. x W 4.5 in. x D 2.75 in.
- Weight: 1 lb. 5 oz.

TORK LCL

Color	Option	Item#
BLK	Adv. w/CG - LE	80-1697
COY	Adv. w/CG - LE	80-1698
OCP	Adv. w/CG - LE	80-1699
ODG	Adv. w/CG - LE	80-1700
RED	Adv. w/CG - LE	80-1701
RGR	Adv. w/CG - LE	80-1702
BLK	Basic w/CG - LE	80-1691
COY	Basic w/CG - LE	80-1692
OCP	Basic w/CG - LE	80-1693
ODG	Basic w/CG - LE	80-1694
RED	Basic w/CG - LE	80-1695
RGR	Basic w/CG - LE	80-1696
BLK	Basic	80-1685
COY	Basic	80-1686
OCP	Basic	80-1687
ODG	Basic	80-1688
RED	Basic	80-1689
RGR	Basic	80-1690

Rx ONLY



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

TORK® LCL (Tactical Operator Response Kit)

REV112624

Basic/Basic with Combat Gauze® Contents:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 6 in. Dressing	 QTY 1 NAR® Gauze, Wound Packing	 QTY 1 NAR® Gauze, Wound Packing	 QTY 1 Combat Gauze® LE
 QTY 1 HyFin® Vent Chest Seal, Twin Pk	 PAIR 2 Responder Black Gloves, LG	 QTY 1 Responder Shears, Lg	-PLUS- BASIC -OR- BASIC with Combat Gauze	

Advanced Contents:

 QTY 1 C-A-T® Tourniquet	 QTY 1 ETD® 6 in. Dressing	 QTY 1 NAR® Gauze, Wound Packing	 QTY 1 Combat Gauze® LE	
 QTY 1 Pre-lubed NAR® Naso, 28F	 QTY 1 HyFin® Vent Chest Seal, Twin Pk	 QTY 1 Enhanced ARS Needle, 10G	 PAIR 2 Responder Black Gloves, LG	 QTY 1 Responder Shears, Lg



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

**Rugged
and Streamlined...**

**Designed to MOLLE to
an assault vest or body armor**