

I-ROK XL *Individual Rugged Outdoor Kit*

REV101024

Essential First-Aid Protected In A Personal Utility Case...

Every good adventure requires that everyone gets home safely. The new I-ROK XL (Individual Rugged Outdoor Kit) from North American Rescue is the perfect companion for any outdoor environment. The I-ROK XL contains the same essential first-aid supplies of the Standard I-ROK Kit with extra supplies to accommodate multiple injuries. These kits are designed to treat common injuries including basic injuries such as cuts, scrapes, bites & stings as well as the life-saving equipment to control major bleeding.

The I-ROK XL comes in the larger ultimate personal R60 utility case that is crushproof, dustproof and adventure proof! Your first-aid and trauma supplies will be protected from the outdoor elements to help ensure your kit is ready if ever needed. Waterproof up to a depth of 2 meters for up to 30 minutes, IP68 rated protection from water, dirt, snow & dust and rubberized protective bumpers to protect against drops. The case handles can be used to secure the kit for your outing as you can strap it to your gear, ATV/UTV, motorcycle, camper or whatever you use to get you to where the adventures are. This ample hard case comes in a variety of colors and keeps your contents protected and secure.

Features:

- Contains essential first-aid supplies to treat basic injuries as well as control major bleeding with extra supplies for multiple injuries
- I-ROK XL features a larger R60 storage case that is crushproof, dustproof and waterproof up to a depth of 2 meters
- Flexible lid organizer helps keep smaller medical supplies neatly stored
- Built-in pressure valve keeps water and dust out while balancing air pressure
- I-ROK XL larger hard case comes in a variety of colors and with the option of Basic, with Bleeding Control Dressing, or with CELOX™ Hemostatic Dressing



Dimensions:

- H 8.9 in. x W 5.4 in. x D 3 in.
- Weight: 3 lb.

I-ROK XL Individual Rugged Outdoor Kit

COLOR	LEVEL	ITEM #
Black	Basic	82-0049
OD Green	Basic	82-0051
Black	with Bleeding Control Dressing	82-0053
OD Green	with Bleeding Control Dressing	82-0055
Black	with Celox™ Hemostatic Dressing	82-0070
OD Green	with Celox™ Hemostatic Dressing	82-0072



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

I-ROK XL Individual Rugged Outdoor Kit

REV101024



Kit Contents:

- 1 x C-A-T® Tourniquet, Orange
- 5 x Antibiotic Ointment
- 1 x NAR® Responder Compressed Gauze
- 3 x Sting & Bite Packet
- 1 x QuikClot® Z-Fold Bleeding Control
- 1 x NAR® Survival Blanket
- 5 x Hand Sanitizer
- 1 x CELOX™ First Aid Hemostatic
- 3 x Salty Britches® Chafing Ointment - Sgl
- 1 x NAR® Mini Responder ETD®, 4in. (Emergency Trauma Dressing)
- 2 pr Responder Nitrile Gloves - Blue, Large
- 1 x HyFin® Vent Compact Chest Seal - Twin Pack
- 1 x NAR® Mini Duct Tape, 2 in.
- 10 x Butterfly Bandage
- 1 x Permanent Marker - Black, Large
- 10 x Bandage, 1 x 3 in.
- 1 x Bleeding Control Instruction Card
- 2 x Gauze Pad, 4 x 4 in. (2 Pk)
- 1 x Triangular Bandage
- 1 x BurnTec® Dressing, 4 x 4 in.
- 1 x Splinter Forceps, 4.5 in.
- 5 x Burn Gel
- 1 x Pelican Ruck R60 Case

The same essential first-aid supplies from our I-ROK Kit with EXTRA supplies to accommodate multiple injuries...



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**