

# NAVAL GUN RESPONSE AID KIT

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



**Designed to meet the rigorous requirements of the U.S. Navy**

For first responder emergency care equipment on board U.S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels. Ergo dynamic, low profile casualty response system features 2 removable equipment pouches that attach via hook-&-loop to a lightweight, durable outer bag, thus meeting the requirements for the U.S. Navy's Naval Gun Response Aid Bag.

Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea. Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/ water resistant, self-repairing zippers with silent 550 cord pulls. Both outer bag & 2 inner pouches are securable using antipilferage seals. Slim-line, lightweight kit stores easily on a wall and can be carried over the shoulder or bandolier-style.



**Pouches easily fit into the standard Naval First Aid Box (not included)**

## FEATURES/SPECIFICATIONS

- Two internal pouches are identical & provide rapid access & organization of critical life saving equipment needed to render on-scene care at sea
- Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self-repairing zippers with silent 550 cord pulls.
- Lightweight, durable slimline pouch is securable using anti-pilferage seals & features a clear window to readily view contents

## NAVAL GUN RESPONSE AID KIT

ITEM #	OPTION	NSN
80-0418	Kit	6545-01-648-2690
80-0417	Bags Only	8465-01-666-7169

## DIMENSIONS

- H 16 in. x W 11 in. x D 5.5 in.
- Weight: 4 lb 11 oz

# NAVAL GUN RESPONSE AID KIT



2 x Internal Pouches


**Each Internal Pouch contains:**



QTY 1 C-A-T® Tourniquet

QTY 1 ETD® 6 in. Dressing

QTY 1 ETD® Abd./Stump



QTY 1 Pre-lubed NAR® Naso, 28F



QTY 1 Naso Airway, 34F



QTY 1 Surgical Lubricant, 5g



QTY 1 CPR Microshield Mask



QTY 1 HyFin® Vent. Twin Pk



QTY 1 Fox Eye Shield w/garter



QTY 1 Surgical Tape, 2 in.



QTY 2 SAM® Splint II



QTY 1 3" Elastic Wrap Bandage



QTY 1 NAR® Trauma Shears, Lg



PAIR 1 Nitrile Gloves, Lg



1 x Naval Gun Response Bag



QTY 1 Permanent Marker, Large

## KIT CONTENTS

- 1x Naval Gun Response Bag, CCRK®
- 2x Internal Pouches (Naval First Aid Box Response Kit)

## INTERNAL POUCH CONTENTS (Each)

- 1x Pouch (CCRK® - Naval First Aid Box Response)
- 1x C-A-T® (Combat Application Tourniquet®)
- 1x ETD® 6 in. Emergency Trauma Dressing
- 1x ETD® Abdominal Emergency Trauma Dressing
- 1x NAR® Nasopharyngeal Airway - 28F (117 MM) - Pre-lubricated
- 1x Nasopharyngeal Airway 34F
- 1x Surgical Lubrication
- 1x CPR Mask, Microshield
- 1x HyFin® Vent Chest Seal Twin Pack
- 1x Fox Eye Shield with Garter
- 1x Roll 2 in. Surgical Tape
- 2x SAM® Splint II
- 1x Elastic Wrap Bandage, 3 in.
- 1x Trauma Shears (7.25 in.)
- 1 pr Bear Claw® Nitrile Gloves, Large
- 1x Permanent Marker, Large

## NAVAL GUN RESPONSE AID KIT

ITEM #	OPTION	NSN
80-0418	Kit	6545-01-648-2690
80-0417	Bags Only	8465-01-666-7169



Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)

