

# MW EVA LITTER 23S

## WITH CASE

### For Particularly Demanding Environments

The MediWear EVA Litter 23S is a quick, easy to assemble lightweight, compact casualty litter that allows rapid transport of casualties. The telescopic poles are constructed out of extremely durable carbon fiber, with new aluminum pole ends that quickly and easily deploy. The litter includes a carrying case that can be secured to load bearing equipment using the MOLLE attachment system.

The light weight and small size of the product make it suitable for carrying into particularly demanding environments where carrying capacity, muscle strength, and reaction time can be the difference between a life lost or a life saved.



### FEATURES/SPECIFICATIONS

- Weighs 6.37 lbs. and can carry up to 330.7 lbs. (165 lbs. per pole)
- Telescopic poles are constructed of carbon fiber with new aluminum carry handles at each pole end, making it significantly easier to deploy
- Has integrated patient belt
- Litter Bed consists of 75% PA and 25% PU
- Fabric finish is water and dirt repellent

### MW EVA LITTER 23S

| ITEM #  | NSN  |
|---------|------|
| 60-0093 | ---- |

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



### DIMENSIONS

- Folded: H 17 in. x W 75 in. x D 5 in.
- Open: H 72.8 in. x W 20 in. x D 2 in.
- Weight of litter (poles and bed): 6.37 lbs.
- Weight (with bag): 7.4 lbs.
- Max carrying load: 330.7 lbs. (165 lbs. per pole)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

# MW EVA LITTER 23S

WITH CASE

**TALON II MODEL 90C COLLAPSIBLE HANDLE LITTER & CARRIER**  
(SOLD SEPERATELY)



**MW EVA LITTER 23S WITH CASE**



FOLDED: H 20.5 in. x W 9 in. x D 10 in.  
OPEN: H 90 in. x W 22.5 in. x D 6 in.  
Weight (with bag): 17 lbs. 13.3 oz.

FOLDED: H 17 in. x W 7.5 in. x D 5 in.  
OPEN: H 72.8 in. x W 20 in. x D 2 in.  
Weight (with bag): 74 lbs.









**1**

**2**

**3**

**4**

**5**

**6**

**7**

*NAR® Heat Reflective Shell™ sold separately*




**NAR-4 Aid Bag sold separately**