

FIELD BLOOD TRANSFUSION KIT - DONOR MODULE KIT

FBTK - DONOR MODULE

“The advantage of the direct transfusion of human blood in cases of severe haemorrhage which we encounter in the emergencies of military surgery cannot to our minds, be overestimated.”¹

The above comment, from a paper presented by Lieutenant Colonel A. Primrose of the Canadian Army Medical Corps in 1916, remains as pertinent today as during those earlier days of military transfusion medicine...

The Field Blood Transfusion Kit (FBTK) from North American Rescue is a comprehensive kit providing all the needed supplies to allow the user to rapidly perform a field blood transfusion. All of the supplies conform with the Joint Theater Trauma System Clinical Practice Guidelines. Equipment provided also allows for verifying the blood types of the donor and recipient as well as the collection of blood samples from the donor for additional surveillance testing.

¹Primrose A, Ryerson ES. **The Direct Transfusion of Blood: Its Value in Haemorrhage and Shock in the Treatment of the Wounded in War.** Br Med J. 1916;2:384–386.

Learn more about this and other **Products with a Mission®** at www.NARescue.com



FEATURES/SPECIFICATIONS

- Comprehensive kit allowing the user to rapidly perform a field blood transfusion
- All supplies are compatible with the Joint Theater Trauma System Clinical Practice Guidelines

- Equipped with Donor/Recipient Modules to verify blood types of the donor and recipient
- Also allows collecting of blood samples from the donor for additional surveillance testing

FBTK - DONOR

ITEM #

83-0001

* Explore our family of blood transfusion products at NARescue.com.

DIMENSIONS

- L 9.5 in x W 14 in
- Weight: 0.75 lbs




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

FIELD BLOOD TRANSFUSION KIT - DONOR MODULE KIT

FBTK - DONOR MODULE

KIT INCLUDES:



KIT CONTENTS

- 1x Donor Placard
- 1x Eclipse Vacutainer Holder w/ Needle, 21ga x 1.25 in.
- 4x Vacutainer, Purple Top (6ml)
- 2x Vacutainer, Red Top (6ml)
- 1x Single Blood Bag (450ml) - DO NOT FOLD
- 1x Blood Bag Sticker
- 1x 550 Cord, 10 in.
- 1x Beaded Cord, 11 in. (Marked at 6.5 in.)
- 1x IV Constricting Band, 1 in. x 18 in.
- 2x Alcohol Prep Pad
- 2x 2-pack Gauze Pad, 4 in. x 4 in.
- 1x Clear Adhesive Dressing, 4 in. x 5 in.
- 2x Antimicrobial Swabsticks, 3 pk
- 1x Kelly Straight Hemostat (Disposable), 4.25 in.
- 1 pr Bear Claw Nitrile Gloves, Lg
- 1x roll of Surgical Tape, 1 in.
- 1x ASBP 572-EWB V29
- 1x Black Permanent Marker

FBTK - DONOR

ITEM #

83-0001

* Explore our family of blood transfusion products at NARescue.com.



Learn more about this and other **Products with a Mission®** at www.NARescue.com

