

# LITTER STANDS (SET OF 2) - 33 IN.

Learn more about this and other **Products with a Mission**<sup>®</sup> at [www.NARescue.com](http://www.NARescue.com)



## Accommodates the Talon<sup>®</sup>, Raven<sup>®</sup>, T3 and Medevac litters

Two, lightweight and compact aluminum folding litter stands, designed to lock into place and provide a stable platform for the patient. Accommodates the Talon, Raven, T3, and Medevac litters to restrict the litter from sliding front to back. Once each stand is opened, the cross member easily locks into place.



## FEATURES/SPECIFICATIONS

- Comes in a set of 2
- Aluminum structure for lightweight strength
- Easily opened to stand 33 in. tall
- Quickly secures the litter to increase stability and access to the patient
- Easily folded for transport and storage
- Accommodates the Talon<sup>®</sup>, Raven<sup>®</sup>, T3 and Medevac litters

## LITTER STANDS (SET OF 2) - 33 IN.

ITEM #	NSN
60-0014	6530-01-523-3314

\* Explore other litter stand options available at [www.NARescue.com](http://www.NARescue.com).

## DIMENSIONS

- Folded: H 34 in. x W 25.75 in. x D 5 in.
- Extended Height: 33 in.
- Base Width: 23.5 in.
- Weight: 8 lb 8 oz
- Unit of Issue: Set of 2
- Payload: 1200 lb



**NORTH AMERICAN RESCUE**<sup>®</sup>  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277