

OSL LITTER STANDS

SET OF 3

Learn more about this and other **Products with a Mission®** at www.NARescue.com

Lightweight and compact aluminum folding litter stands

Designed to lock into place and provide a stable platform for the patient, the OSL litter stands come in a set of 3 and accommodate the oversized OSL litter. Equipped with litter stirrups to restrict the litter from sliding front to back. Once the stands are opened, the cross members easily lock into place.



FEATURES/SPECIFICATIONS

- Set of 3 lightweight aluminum stanchions to accommodate our oversized OSL™ Litter
- Added stability for the patient
- Better access for the medic
- Easy to open; quickly locks to secure litter
- Easily folded for transport & storage

OSL LITTER STANDS SET OF 3

ITEM #	NSN
60-0029	----

DIMENSIONS

- Folded: H 34 in. x W 25 in. x D 2.5 in.
- Extended Height: 32 in.
- Base Width: 23.5 in.
- Weight: 28.5 lb
- Unit of Issue: Set of 3
- Payload: 1600 lb