

QUIKLOT® EMS 4 X 4

HEMOSTATIC GAUZE

Learn more about this and other **Products with a Mission®** at www.NARescue.com



Effective, simple and safe

QuikClot® EMS 4 x 4 dressing is a soft, white, sterile, non woven 4 in. by 4 in. 4-ply gauze impregnated with the hemostatic, kaolin, an inert mineral that does not generate heat. Each QuikClot® EMS 4 X 4 dressing is individually wrapped in an easy rip, pouch. Indicated for temporary external control of traumatic bleeding, QuikClot® EMS 4 X 4 dressing is flexible and pliable and contours to all wounds. Effective, simple and safe to use when trying to control life-threatening hemorrhage from wounds that are not amenable to tourniquet application.



FEATURES/SPECIFICATIONS

- **Controls Bleeding:** Promotes clotting within minutes, allowing responders to control bleeding in an efficient and timely manner.
- **Easy To Use:** QuikClot® EMS 4 x 4 Hemostatic Gauze uses the same placement technique as standard gauze
- **Heat Free:** QuikClot® products produce NO Heat

QUIKLOT® EMS - 4X4

ITEM #

30-0095

DIMENSIONS

- Packaged - H 5.5 in. x W 5 in. x D 0.2 in.
- Gauze - W 4 in. x L 4 in.
- Weight: 0.8 oz

* Search other QuikClot® products at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277