

TCCC MAN SIMULATOR

Learn more about this and other **Products with a Mission®** at www.NARescue.com



Enhances TCCC Instruction

Rugged, state-of-the-art TCCC MAN® training simulator provides enhanced training in Tactical Combat Casualty Care (TCCC) from Point of Injury to transfer of higher care!

ALL-IN-ONE, hyper-realistic trauma simulation solution enabling instructors to build multiple realistic traumatic casualty scenarios with reactive patients ideal for use in classroom settings, trauma lanes and field training scenarios. Specifically **designed for the application of skills and education** utilizing protocols from **Tactical Combat Casualty Care (TCCC)**



Includes GSW chest with anterior entry and mid-axillary exit wound, deep packable pelvic wound, complete below-the-knee amputation let (TQ responsive), and unstable pelvis with palpatus crepitus



FEATURES/SPECIFICATIONS

- Remote controlled bleeding from multiple sites
- Arterial and venous bleeding
- Quick-exchange Lithium-ion battery
- Large capacity on-board bleeding bag
- Weight 150 lbs
- Industry-leading durability -withstands drops, falls, and drags
- Welded steel frame, steel joint and bones
- Austere environmental use (rain, snow, mud, sand)
- Available in Light Tone or Deep Tone skin
- Includes 3-year manufacturer warranty from date of delivery

TCCC MAN®

ITEM #	NSN
93-0958	----

* Explore more simulator options at www.NARescue.com

DIMENSIONS

- L 40 in. x H 36 in. x W 48 in.
- Weight: 150 lbs



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

TCCC MAN SIMULATOR

KIT CONTENTS

- 1x Remote instructor tablet
- 1x Blood reservoir pressure fill tank
- 1x Velcro digital pattern uniform
- 2x 3-wheeled deployment-style bags for safe storage and transport
- 1x NAR Training IFAK
- 1x Operator's Manual

TCCC INTERVENTIONS SUPPORTED

- Massive hemorrhage control (tourniquet application and wound packing)
- Chest seal application
- Pleural needle decompression (multiple sites)
- Airway management (NP/OP airways, supraglottic devices, endotracheal intubation)
- Tibial intraosseous site

INJURIES INCLUDED

- GSW chest with anterior entry and mid-axillary exit wound
- Deep packable pelvic wound
- Complete below the knee amputation leg (TQ responsive)
- Unstable pelvis with palpable crepitus

ADD'L SPECS

- Features proprietary communication link that requires NO Wi-Fi, Bluetooth, or LAN connections
- Simultaneous long-range communications with numerous NARS devices



TCCC MAN®

ITEM #	NSN
93-0958	----

Learn more about this and other Products with a Mission® at www.NARescue.com

