

CHEST TUBE INSERTION KIT

The NAR® Chest Tube Insertion Kit is designed to provide the necessary equipment needed to perform emergency chest tube thoracostomy. The compact and lightweight kit is packaged in an easy to open, vacuum sealed pouch, the kit can be stored in your kit or crash cart.



Learn more about this and other Products with a Mission® at www.NARescue.com



all KITS INCLUDE:

 QTY 1 Petrolatum Gauze	 QTY 1 Protected #10 Scalpel	 QTY 1 5.5" Kelly Hemostat, Cvd	 QTY 1 Surgical Latex Gloves	 QTY 1 Chest Tube Valve
 QTY 1 Chloraprep Swabsticks, 3Pk	 QTY 1 4" x 4" Gauze Pad 2 PK	 QTY 1 O-Silk Sutures, 30 in.	 QTY 1 Roch-Peon, Cvd - Sterile, 8.5"	

+plus, by kit type:

28FR	 QTY 1 Chest Tube 28 French, Soft	- Or -	36FR	 QTY 1 Chest Tube 36 French, Soft
-------------	---	--------	-------------	---

FEATURES/SPECIFICATIONS

- Designed to provide the necessary equipment needed to perform emergency chest tube thoracostomy
- Compact and lightweight kit packaged in an easy to open, vacuum sealed pouch
- Can be stored in your kit or crash cart

KIT CONTENTS

- 1 x Petrolatum Gauze - 3 in. x 18 in.
- 1 x Protected #10 Scalpel
- 1 x Hemostat, Kelly Curved - 5.5 in.
- 1pr Sterile Surgical Latex Gloves - 8.5 in.
- 1 x Chest Tube Valve
- 1 x Chest Tube, 28 **or** 36 French Soft
- 1 x Chloraprep Swabsticks, 3 per pkg
- 1 x Gauze Pad (2 pk) - 4 in. x 4 in.
- 1 x O-Silk Suture, 30 in.
- 1 x Curved Roch-Peon, Sterile - 8.5 in.

CHEST TUBE INSERTION KIT

ITEM #	SIZE	NSN
83-0007	36FR	6515-01-706-5310
83-0039	28FR	----

** There are other options available at a place for this item and other items.*

DIMENSIONS

- L 12 in. x W 7.5 in.
- Weight: 9.8 oz