

AMPHIBIOUS AID BASIC KITS

WATERTIGHT

Keep Sensitive Gear Dry While Conducting Maritime And Underwater Operations

The Amphibious Aid Kit (Watertight)- from North American Rescue® is designed as a portable First Aid Kit for use in the workplace that exceeds the ANSI/ISEA Z308.1-2021 Standard for Class B kits. This is accomplished by adding additional bleeding control elements proven to save lives. To bring this kit to the market, NAR® partnered with Watershed®, the manufacturer of drybags and other watertight containers used by the US military special operations community to keep their sensitive gear dry while conducting maritime and underwater operations. The kit contents are arranged in labeled and zipper-accessed mesh pouches. This kit supplies not only the item needed to treat the most common workplace injuries but many of the most serious, life threatening with the C-A-T® Tourniquet, the safest and most effective prehospital tourniquet on the market today. This kit comes in Basic, Basic with Combat Gauze®, and Basic with Celox™.



Learn more about this and other Products with a Mission® at www.NARescue.com



AMPHIBIOUS AID BASIC KITS - (WATERTIGHT)

ITEM #	COLOR	OPTION	NSN
80-1559	Orange	Basic	----
80-1560	Ranger Green	Basic	----
80-1561	Orange	Basic with Combat Gauze®	----
80-1562	Ranger Green	Basic with Combat Gauze®	----
80-1621	Orange	Basic with Celox™ Gauze	----
80-1622	Ranger Green	Basic with Celox™ Gauze	----

* Explore advanced versions of this kit at www.NARescue.com.

FEATURES/SPECIFICATIONS

- A portable First Aid Kit for use in the workplace that exceeds the ANSI/ISEA Z308.1-2021 Standard for Class B kits.
- NAR® partnered with Watershed to keep their sensitive gear dry while conducting maritime and underwater operations.
- Supplies not only the item needed to treat the most common workplace injuries but many of the most serious, life threatening ones
- Comes in Basic, Basic with Combat Gauze®, and Basic with Celox™

DIMENSIONS

- H 6.5 in. x W 9 in. x L 11.5 in.
- Weight: 4 lbs 10 oz


AMPHIBIOUS AID BASIC KITS

WATERTIGHT

all KITS INCLUDE:



+plus, by kit type:



KIT CONTENTS

- 2 x C-A-T® - Orange (Combat Application Tourniquet®)
 - 2 x ETD® 4 in Flat Responder Emergency Trauma Dressing
 - 1 x HyFin® Vent Compact Chest Seal Twin Pack
 - 2 x NAR® Wound Packing Gauze - Z-Fold 3 in x 5 yd Dressing
 - 2 x Elastic Wrap Bandage - 4 in 5 yd
 - 1 x Triangular Bandage - 40 in x 40 in x 56 in
 - 2 x Pad Gauze - 4 in x 4 in (2 Pack)
 - 2 x Pad Gauze - 5 in x 9 in
 - 1 x Instant Cold Pack - 6 in x 9 in
 - 1 x Universal Malleable Splint
 - 1 x NAR® Survival Blanket - 56 in. x 96 in.
 - 5pr Sand Bear Claw® Nitrile Glove Kit
 - 1 x 7.25 in. Responder Trauma Shears
 - 2 x Sterile Oval Eye Pad - 1.625 in x 2.625 in
 - 1 x Gecko Grip Multi-purpose Flat Tape (6 Count Strips - 2" x 15")
 - 1 x Elastic Adhesive Tape - 4 in 5 yd
 - 1 x Large Permanent Marker, Black
 - 1 x MED Cross Patch (2"x2") - Black/Green
- +plus**
- 2 x NAR® Wound Packing Gauze - Z-Fold 3 in x 5 yd Dressing
 - or -
 - 2 x Combat Gauze® - Z-fold 3 in x 4 yd Dressing
 - or -
 - 2 x Celox™ Gauze - z-fold 3 in x 5ft Dressing



Learn more about this and other Products with a Mission® at www.NARescue.com



AMPHIBIOUS AID BASIC KITS - (WATERTIGHT)

ITEM #	COLOR	OPTION	NSN
80-1559	Orange	Basic	----
80-1560	Ranger Green	Basic	----
80-1561	Orange	Basic with Combat Gauze®	----
80-1562	Ranger Green	Basic with Combat Gauze®	----
80-1621	Orange	Basic with Celox™ Gauze	----
80-1622	Ranger Green	Basic with Celox™ Gauze	----

* Explore advanced versions of this kit at www.NARescue.com.