

# HEADREST IFAK KIT

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



## Optional Action Statement or Major Fact

The Headrest IFAK Kit was designed to provide personnel with an easily accessible and quick to deploy medical pouch that mounts easily to a vehicle's headrest. The large elastic band keeps the kit stowed securely, while the adjustable attachment straps keep the kit tight against a headrest. The durable side pull handles allow for the kit bag to be rapidly deployed from either side of the mount. The 550 cord zipper pulls are reliable and allow for quick access to the bag's contents.

This Kit comes fully loaded with critical medical equipment for point-of-wounding treatment during traumatic injuries. The design of the Headrest IFAK platform features durable nylon construction, multiple interior elastic loops allowing for the contents to be organized and secured.



## FEATURES/SPECIFICATIONS

- Conveniently mounts to most vehicle headrests
- Elastic band keeps the kit bag stowed securely
- Easy to grab side handles allow for easy deployment
- IFAK pouch features durable 550 cord pulls

- The pouch design allows for the contents of an individual first aid kit or vacuum sealed kit
- Multiple interior elastic loops for secure gear organization
- Provide emergency treatment solutions until definitive care is available

## HEADREST IFAK KIT

ITEM #	COLOR	OPTION	NSN
80-1091	Black	Basic	----
80-1150	Coyote	Basic	----
80-1151	Black	Basic with Combat Gauze®	----
80-1152	Coyote	Basic with Combat Gauze®	----

## DIMENSIONS

- H 5.25 in. x W 7.25 in. x D 4.25 in.
- Weight: 1.62 lbs




# HEADREST IFAK KIT

**all KITS INCLUDE:**



**+plus, by kit type:**



## KIT CONTENTS

- 2 x Combat Application Tourniquet, (C-A-T®) - Orange
- 1 x NAR® Wound Packing Gauze - 3 in. x 5 yds.
- 2 x NAR® Mini Responder ETD®, 4 in. (Emergency Trauma Dressing)
- 2 x HyFin® Vent Chest Seal - Twin Pack
- 1 x Triangular Bandage - 40 in. x 40 in. x 56 in.
- 20 x Bandage - 1 in. x 3 in.
- 2 pr Black Responder glove kit - Large
- 1 x Microshield CPR mask
- 1 x NAR® Survival Blanket

**+plus**

- 1 x NAR® Wound Packing Gauze - 3 in. x 5 yds. - *or* -
- 2 x Combat Gauze® Z-fold Hemostatic Dressing



## HEADREST IFAK KIT

ITEM #	COLOR	OPTION	NSN
80-1091	Black	Basic	----
80-1150	Coyote	Basic	----
80-1151	Black	Basic with Combat Gauze®	----
80-1152	Coyote	Basic with Combat Gauze®	----

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

