

# X-PAK

## Expeditionary Personnel Aid Kit

Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



### Address the first aid needs of a small group

From the mountains to the beach—no adventure should be without the X-PAK. NAR's Expeditionary Personnel Aid Kit (X-PAK) equips you to address the first aid needs of a small group with injuries ranging from major bleeding trauma and fractures to minor cuts and abrasions. The compact size enables optimized stowage for easy transport, ensuring that you can always have the necessary equipment to save a life within reach.

#### Ideal For

- Individual Sports
- Backpacking
- Hiking
- Hunting
- Camping
- Vehicles
- Adventure travel
- Surfing
- Motocross / ATV
- Mountain Biking
- Skiing
- Snowboarding



**First Aid for minor to severe injuries for multiple patients to provide emergency treatment until definitive care arrives**

### FEATURES/SPECIFICATIONS

- Treat minor to severe injuries – 2 to 3 patients
- Provide emergency treatment solutions until definitive care is available
- Compact and capable – ideal for adventure activities or extended travel

### X-PAK

ITEM #	NSN
82-0001	----

### DIMENSIONS

- L 9 in x W 5 in x H 6.5 in
- Weight: 2.72 lb




**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# X-PAK

## Expeditionary Personnel Aid Kit

**all KITS INCLUDE:**



### KIT CONTENTS

#### MAJOR TRAUMA

- 2 x C-A-T® - ORG (Combat Application Tourniquet®)
- 1 x 6 in. Responder ETD® (Emergency Trauma Dressing)
- 1 x Compressed Gauze (4.5 in. x 4.1 yd)
- 1 x QuikClot® Bleeding Control Dressing
- 2 x NAR® Survival Blanket

#### BANDAGES

- 2 x Burntec® Burn Dressing, 4 in x 4 in (10 cm x 10 cm)
- 20 x Flexible Fabric Bandage, 0.75 in x 3 in
- 2 x Wound Closure Strip (Pack of 10)
- 2 x Gauze Pad, 4 in x 4 in
- 2 x Clear Adhesive Dressing, 4 in x 5 in
- 2 x Clear Adhesive Dressing, 2.4 in x 2.75 in
- 1 x Silver Bandage, 2 in x 2 in
- 1 x Silver Bandage, 4 in x 4.5 in
- 2 x Petrolatum Gauze, 3 in x 18 in
- 1 x NAR® Flat Tape

#### WOUND CARE

- 1 x Irrigation Syringe, 10cc
- 1 x Saljet Rinse Cleanser

#### SPLINTING

- 1 x Triangular Bandage, 40 in x 40 in x 56 in

#### EYE CARE

- 2 x Polycarbonate Eye Shield

#### PERSONAL PROTECTION AND ACCESSORIES

- 3 pr Bear Claw Nitrile Trauma Gloves-LG
- 1 x Trauma Shears (LG), 7.25 in
- 1 x Permanent Marker- Blk, Small
- 1 x Penlight
- 1 x CPR Microshield Mask
- 1 x Splinter Forceps, 4.5 in
- 1 x Luminous MED ID Patch



### X-PAK

ITEM #	NSN
82-0001	----

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

