

CRISIS INCIDENT RESPONSE KITS

Learn more about this and other **Products with a Mission**® at www.NARescue.com



Treat Several Injured Persons

The Crisis Incident Response Kit from North American Rescue is designed to provide the necessary medical products to treat the more serious medical emergencies that may occur in a school setting. This ruggedized Nylon kit was designed to combat life threatening emergencies such as inability to breathe or major bleeding. Available in black and high-visibility orange, this kit provides carrying options such as over the shoulder, hand carry and fanny-pack style.

FEATURES/SPECIFICATIONS

- The necessary emergency first aid products to treat several injured persons ranging from minor to severe injuries
- The School First Aid Kit designed to treat penetrating trauma and major bleeding as well as the administration of rescue breathing & basic first aid
- High visibility orange bag constructed of durable 1000 denier nylon with reflective trim
- Multiple carry options including over the shoulder, hand carry & fanny pack style

CRISIS INCIDENT RESPONSE KITS

ITEM #	COLOR
85-0410	Orange
85-1338	Black

DIMENSIONS

- Packaged: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 5 lb

CRISIS INCIDENT RESPONSE KITS

all KITS INCLUDE:



KIT CONTENTS

- 2 x C-A-T® Tourniquet - Orange or Black
- 4 x NAR® Compressed Gauze
- 4 x NAR® 4 in. Flat Emergency Trauma Dressing (ETD®)
- 4 x NAR® 6 in. Emergency Trauma Dressing (ETD®)
- 1 x NAR® Abdominal Emergency Trauma Dressing (ETD®)
- 2 x HyFin® Vent Compact, Twin Pack
- 2 x CPR Microshield
- 2 x SAM® Splint II
- 1 x roll 2 in. Surgical Tape
- 5 pr Bear Claw Gloves, large
- 1 x NAR® 7.25 in. Trauma Shears
- 50 x Flexible Fabric Bandages
- 2 x NAR® Survival Blanket - 56 in. x 96 in.



Learn more about this and other Products with a Mission® at www.NARescue.com



CRISIS INCIDENT RESPONSE KITS

ITEM #	COLOR
85-0410	Orange
85-1338	Black